



## Bertec® Portable Essential

The Bertec® Portable Essential™ provides a variety of assessments, including standard static balance tests, weight shifting, and optional higher-level balance testing for athletes (COBALT™). Static balance and weight shifting activities can be easily integrated into a customized therapy program.

- High-sensitivity dual-balance force plate
- Objective assessment of balance with normative data for comparison
- Convenient for clinical rehabilitation, home health, bedside, and community health programs
- Portable to bring balance services beyond clinic walls
- Training menus for use with biofeedback rehabilitation
- Low profile and light weight



**Standard Package**

Bertec Balance Advantage Software	Foam balance pad (19" x 20" x 4")
Dual-balance force plate	Database merge and sync
Laptop computer with touchscreen display and remote control	One-year parts and labor with ongoing support
Power cable and USB cable	
Wireless pointing device/mouse	

**Assessment Protocols**

Modified Clinical Test of Sensory Interaction on Balance (mCTSIB)

Limits of Stability (LOS)

Rhythmic Weight Shift (RWS)

Unilateral Stance (US)

Weight Bearing Squat (WBS)

**Featured Training**

Quick Training

Mobility

Closed Chain

Weight Shifting

Seated Training

**Specifications**

High-sensitivity dual-balance force plate

Standardized protocols with normative data

Low 1.5" profile in an 18" x 20" configuration

Modular software design to fit your needs

500 lb load capacity

Patented and HIPAA compliant patient database with merge and sync capabilities across all Bertec® Balance Advantage® products

Portable system powered by USB

**Options**

Bertec® Vision Advantage™

COBALT™ powered by Bertec®

Carrying case (custom fit to all components)

Extended warranty

For more information, contact Bertec at 614-543-8099 or by email at [info@bertec.com](mailto:info@bertec.com)