



## Pressure Prints


Average



Maximum



## COP Diagram



## COP Parameters

Length of gait line, mm	Left	76±12	
	Right	76±11	
	Diff, %	0.2	
Single support line, mm	Left	76±12	
	Right	76±11	
	Diff, %	-0.1	
Ant/Post position, mm		178±3	
Lateral symmetry, mm		3±6	




**Patient**  
 Project  
 First Name  
 Last Name  
 Sex

Project 1  
 Demo Records Treadmill Re...  
 MyoPressure-T  
 Male

**Record**  
 Name  
 Date Measured  
 Number of periods

Subject 12mph bf - Runnin...  
 20.08.2013 02:49  
 39



## Running Phase Parameters

Stance phase, %	Left	26.2±0.6	
	Right	25.3±0.7	
	Diff, %	-3.3	
Swing phase, %	Left	73.8±0.6	
	Right	74.7±0.7	
	Diff, %	1.2	
Initial swing, %	Left	24.1±0.8	
	Right	24.5±0.9	
	Diff, %	1.7	
Mid swing, %	Left	25.3±0.8	
	Right	26.2±0.6	
	Diff, %	3.3	
Terminal swing, %	Left	24.5±0.8	
	Right	24.1±0.7	
	Diff, %	-1.7	
Stance duration, ms	Left	162±4	
	Right	157±5	
	Diff, %	-3.2	
Swing duration, ms	Left	456±7	
	Right	462±8	
	Diff, %	1.2	
Total flight duration, ms		300±8	

## Running Spatial Parameters

Foot rotation, deg	Left	6.2±2.0	
	Right	6.8±1.3	
Step length, cm	Left	165±2	
	Right	162±2	
	Diff, %	-1.7	
Stride length, cm		327±3	
Step width, cm		4±2	
Velocity, km/h		19.0±0.2	

## Running Time Parameters

Step time, ms	Left	308±7	
	Right	311±7	
	Diff, %	0.9	
Stride time, ms		618±7	
Cadence, step/min		194±2	




**Patient**  
 Project Project 1  
 First Name Demo Records Treadmill Re...  
 Last Name MyoPressure-T  
 Sex Male


**Record**  
 Name Subject 12mph bf - Runnin...  
 Date Measured 20.08.2013 02:49  
 Number of periods 39




## Left-Right Curve Comparison




## Minimum, N




## Maximum, N



Marker, 9.86 s



### Pressure distribution, Gait




### Maximum Force, N

Force LT	1237.68	
Force RT	1179.03	


### Mean Force, N

Force LT	0.00	
Force RT	522.98	

Marker, 9.94 s



### Pressure distribution, Gait



### Maximum Force, N

Force LT	1237.68	
Force RT	1179.03	


### Mean Force, N

Force LT	0.00	
Force RT	1170.27	

Marker, 10.19 s



### Pressure distribution, Gait




### Maximum Force, N

Force LT	1135.41	
Force RT	1219.09	


### Mean Force, N

Force LT	816.47	
Force RT	0.00	

Marker, 10.25 s



### Pressure distribution, Gait



### Maximum Force, N

Force LT	1135.41	
Force RT	1219.09	

### Mean Force, N

Force LT	1135.41	
Force RT	0.00	




**Patient**  
Project Project 1  
First Name Demo Records Treadmill Re...  
Last Name MyoPressure-T  
Sex Male

**Record**  
Name Subject 12mph bf - Runnin...  
Date Measured 20.08.2013 02:49  
Number of periods 39



Marker, Statistics

Lenovo EasyCamera



Pressure distribution, Gait

